



Hapo Exoskeleton Cleaning Guide

To ensure the continued performance and hygiene of your Hapo Exoskeleton, regular cleaning is essential. Please follow the guidelines below for proper care and maintenance.

Routine cleaning: after each use

The Hapo Exoskeleton should be wiped down after each use to maintain its cleanliness. Use a **quat-based disinfectant spray** (quaternary ammonium) that does not contain hydrogen peroxide, bleach, or alcohol. (e.g., Surfa-safe premium). These disinfectants are safe for all materials used in the exoskeleton and effectively disinfect without damaging the product.

- **Recommended use:** Spray the disinfectant on a soft cloth and wipe all exposed surfaces or apply directly to the device.
- **Frequency:** This cleaning is sufficient for routine maintenance after each use as long as there are **no visible signs of soiling or contamination**.

Extensive cleaning: cleaning with Visible Dirt or Foreign Material

In the event that visible dirt or foreign material is present on the exoskeleton, further cleaning will be required. Follow these steps:

1. **Disassemble** the exoskeleton as per the instructions in the user manual.
2. **Wash the textiles:** The fabric components should be washed according to the guidelines provided in the user manual. This will ensure the integrity and cleanliness of the fabric.
3. **Wipe down hard materials:** Use a clean cloth dampened with water and pH-neutral soap to gently wipe down the hard components of the exoskeleton.
4. **Wipe dry immediately:** Ensure all parts are thoroughly dried before reassembling the exoskeleton to prevent any moisture build-up.

Extensive Cleaning should take place at least once a month to ensure the exoskeleton is maintained in best conditions.

By following these simple cleaning steps, you will ensure that your Hapo Exoskeleton remains in optimal condition for everyday use.

For more detailed instructions on disassembly and deep cleaning, please refer to the **User Manual**. Regular maintenance is crucial for ensuring the safety, comfort, and longevity of your exoskeleton.